




FRONTAL, SUPERIOR AND POSTERIOR VIEWS SHOWING STANDARD ELECTRODE POSITIONS

T3/T4: 1/2" above ears

CZ: Intersection of midlines

Fp1/Fp2: 1" from midline

O1/O2: 1" from midline

C3: half way between T3 & CZ

C4: Halfway between T4 & CZ

F3: Halfway between C3 & Fp1

F4: Halfway between C4 and Fp2

P3: Halfway between C3 & O1


P4: Halfway between C4 & O2



Equals 1" from Midline



Horizontal Circumference around head. Starts and ends at the Temporal sites



SINGLE PLANE PROJECTION OF HEAD SHOWING STANDARD ELECTRODE POSITIONS

Figure 14. 10-20 Electrode Placement System.

10-20 Treatment Locations